



124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536-4445

888.395.2221 / 508.540.3444
custservice@acciusa.com

26 Aug 2026

Subject: Gen5 Protocol for Pyros eXpress® Software User Validation Protocol with Agilent BioTek Epoch 2 Absorbance Microplate Spectrophotometer

This notification is to inform you of updates to the Gen5 protocol file provided to execute Test case 6.5 in Section 6 of the Epoch 2 User Validation Protocol (G_3045). Protocol file TC_6.5 rev3 is uploaded to ACC Digital Distribution Portal (Associates of Cape Cod, Inc.).

The updates to TC_6.5 rev3 are as follows:

- This Gen5 protocol file now requires Gen5 version 3.18 or higher and is not reverse compatible with lower versions of Gen5 (Agilent software releases can be downloaded at no cost from Agilent Software Download at [BioTek Software Releases | Agilent](#)).
- This Gen5 protocol file is locked to prevent unauthorized changes.
- The data reduction was updated in this protocol file to more closely align with Agilent specifications for evaluation of a water baseline check. There is no impact on data previously evaluated using previous revisions of this protocol.

We appreciate your continued partnership. Please do not hesitate to contact us at techservice@acciusa.com with any questions.

Sincerely,
Sean A. Hauber
Marketing Product Manager
SHAuber@acciusa.com
858-356-7600

MKT26-353